

Pinocembrin

Catalog No: #S2992


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pinocembrin
Brief Description	Inhibitors
Purification	99.50%
Target Name	NF-κB inhibitor
Calculated MW	256.25
Formulation	C15H12O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
 SMILES: C1[C@H](OC2=CC(=CC(=C2C1=O)O)O)C3=CC=CC=C3
 Pathways: NF-κB Receptor: NF-κB
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Saad MA, et al. Pharmacol Rep. 2015 Feb;67(1):115-22.

Note: This product is for in vitro research use only