

Cordycepin

Catalog No: #S2993


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cordycepin
Brief Description	Inhibitors
Purification	99.50%
Target Name	TNF inhibitor
Calculated MW	251.24
Formulation	C10H13N5O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: C1[C@H](O[C@H]([C@@H]1O)N2C=NC3=C2N=CN=C3N)COP
Pathways: Apoptosis
Receptor: TNF- α
Boiling pt: Melting pt: Solubility: DMSO: Soluble
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Yang J, et al. Int J Immunopathol Pharmacol. 2015 Sep;28(3):296-307.

Note: This product is for in vitro research use only