

## Dihydromyricetin

Catalog No: #S2998


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Dihydromyricetin
Brief Description	Inhibitors
Purification	96.68%
Target Name	GABA Receptor agonist
Calculated MW	320.25
Formulation	C15H12O8
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: C1=C(C=C(C(=C1O)O)O)[C@@H]2[C@H](C(=O)C3=C(C=C(C=C3O2)O)O)O  
 Pathways: Membrane transporter/Ion channel; Neuroscience  
 Receptor: GABAR  
 Boiling pt: Melting pt: Solubility: DMSO: 64mg/ml(199.84 mM); Water:

## References

1. Zhang YS, et al. Yao Xue Xue Bao, 2003, 38(4), 241-244.

Note: This product is for in vitro research use only