

## Batyl alcohol


Catalog No: #S3004

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Batyl alcohol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	344.58
Formulation	C21H44O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:CCCCCCCCCCCCCCCCOCC(O)COPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Suemune H, et al. Chem Pharm Bull (Tokyo). 1987 Aug;35(8):3112-8.

Note: This product is for in vitro research use only