

4-Hydroxychalcone

Catalog No: #S3005


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	4-Hydroxychalcone
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	224.26
Formulation	C15H12O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=CC=C(C=C1)C(=O)/C=C/C2=CC=C(C=C2)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Qu Q, et al. Biomed Res Int. 2014;2014:603415.

Note: This product is for in vitro research use only