

5-Phenyl-2,4-pentadienoic acid

Catalog No: #S3006


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	5-Phenyl-2,4-pentadienoic acid
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	174.2
Formulation	C11H10O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: C1=CC=C(C=C1)/C=C/C=C/C(=O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Werbel LM, et al. J Med Chem. 1968 Sep;11(5):1073-4.

Note: This product is for in vitro research use only