

Spermine

Catalog No: #S3007


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Spermine
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	202.35
Formulation	C10H26N4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:NCCCCCCCCNCCCNPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Elustondo PA, et al. Biochim Biophys Acta. 2015 Feb;1847(2):231-40.

Note: This product is for in vitro research use only