

2-Methoxy-1,4-naphthoquinone


Catalog No: #S3008

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | 2-Methoxy-1,4-naphthoquinone |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | JNK; p38 MAPK inhibitor |
| Calculated MW | 188.18 |
| Formulation | C11H8O3 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Immunology
SMILES: COC1=CC(=O)C2=CC=CC=C2C1=O
Pathways: MAPK Receptor: JNK; p38 MAPK
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ong JY, et al. Life Sci. 2015 Aug 15;135:158-64.

Note: This product is for in vitro research use only