

Bovinoicidin


Catalog No: #S3009

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|-------------------------------------------------|
| Product Name | Bovinoicidin |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | Others |
| Calculated MW | 119.08 |
| Formulation | C3H5NO4 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
SMILES: OC(=O)CC[N+](=O)[O-]
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Gupta S, Sharma B. Eur J Pharmacol. 2014 Jun 5;732:111-22.

Note: This product is for in vitro research use only