

Xanthurenic Acid

Catalog No: #S3010


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Xanthurenic Acid
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	205.17
Formulation	C10H7NO4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: OC(=O)C1=CC(=O)C2=C(N1)C(O)=CC=C2
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Haruki H, et al. J Biol Chem. 2015 Nov 12. pii: jbc.C115.680488.

Note: This product is for in vitro research use only