

## Kynurenine


Catalog No: #S3011

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kynurenine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	208.22
Formulation	C10H12N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: NC(CC(=O)C1=CC=CC=C1N)C(=O)O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: Yellow Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Liu JJ, et al. *Biochim Biophys Acta*. 2015 May;1852(5):980-91.

Note: This product is for in vitro research use only