

Mangiferin

Catalog No: #S3012


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Mangiferin |
| Brief Description | Inhibitors |
| Purification | 95.00% |
| Target Name | Others |
| Calculated MW | 422.35 |
| Formulation | C19H18O11 |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research

Area: Other
 SMILES: C1=C2C(=CC(=C1O)O)OC3=CC(=C(C(=C3C2=O)O)[C@H]4[C@@H]([C@H]([C@@H]([C@H](O4)CO)O)O)O)O
 Pathways: Other
 sReceptor: Others
 Boling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: Yellow Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Miura T, et al. Phytomedicine. 2001 Mar;8 (2):85-7.

Note: This product is for in vitro research use only