

## Catharanthine tartrate

Catalog No: #S3013


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Catharanthine tartrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	486.51
Formulation	C25H30N2O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CCC1=C[C@H]2C[C@]3([C@@H]1N(C2)CCC4=C3NC5=CC=CC=C45)C(=O)OC.[C@@H]([C@H](C(=O)O)O)(C(=O)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Lin C, et al. Biomed Chromatogr. 2015 Jan;29(1):97-102.

Note: This product is for in vitro research use only