

Rheochrysidin

Catalog No: #S3019


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rheochrysidin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	284.27
Formulation	C16H12O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:COC1=CC(O)=C2C(=O)C3=C(O)C=C(C)C=C3C(=O)C2=C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Ethanol: 1 mg/mL (3.51 mM); DMSO: 6 mg/mL (21.1 mM)Appearance:Orange PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Ma X, et al. Pest Manag Sci, 2010, 66(7), 718-724.

Note: This product is for in vitro research use only