

## Sclareol

Catalog No: #S3020


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sclareol
Brief Description	Inhibitors
Purification	99.58%
Target Name	Others inhibitor
Calculated MW	308.51
Formulation	C <sub>20</sub> H <sub>36</sub> O <sub>2</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CC(C)(C)C1CCC(C)(C)C2C(C)C(C)C(C)C12  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: Ethanol: 24 mg/mL warmed (77.79 mM); DMSO: 61 mg/mL warmed (197.73 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Dimas K, et al. Biomed Pharmacol, 2006, 60(3), 127-133.

Note: This product is for in vitro research use only