

## Bulleyaconi cine A

Catalog No: #S3021


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Bulleyaconi cine A
Brief Description	Inhibitors
Purification	98.00%
Target Name	Sodium Channel inhibitor
Calculated MW	627.76
Formulation	C35H49NO9
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: CCN1C[C@@]2(CC[C@@H]([C@@]34[C@@H]2[C@H]([C@H](C31)[C@]5(C[C@@H]([C@]6(C[C@@H]4[C@@H]5[C@H]6C(=O)C7=CC=C(C=C7)OC)O)OC(=O)C)OC)COC Pathways: Membrane transporter/Ion channel Receptor: Sodium Channel Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Solid Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Wang CF, et al. Anesth Analg. 2008 Oct;107(4):1397-405.

Note: This product is for in vitro research use only