

## Sparteine sulfate

Catalog No: #S3029


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sparteine sulfate
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others
Calculated MW	332.46
Formulation	C15H28N2O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:OtherSMILES:OS(O)(=O)=O.C1CCN2CC3CC(CN4CCCC34)C2C1Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Silva MR, et al. Genet Mol Res. 2014 Dec 12;13(4):10510-7.

Note: This product is for in vitro research use only