

L-Sophoridine

Catalog No: #S3030


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	L-Sophoridine
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	248.36
Formulation	C15H24N2O
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:O=C1CCCC2C3CCCN4CCCC(CN12)C34Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:10 mM in DMSOAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Liu J, et al. 2015 May;31(5):585-9, 595. Chinese.

Note: This product is for in vitro research use only