

Deoxycorticosterone acetate

Catalog No: #S3033


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Deoxycorticosterone acetate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	372.5
Formulation	C23H32O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CC(=O)OCC(=O)[C@H]1CC[C@@H]2[C@@]1(CC[C@H]3[C@H]2CC=C4[C@@]3(CC[C@@H](C4)O)C)C Pathways: Others Receptor: Others
 Boiling pt: Melting pt: Solubility: Ethanol: 6 mg/mL (16.1 mM); DMSO: 11 mg/mL (29.53 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Schenk J, et al. J Pharmacol Toxicol Methods, 1992, 27(3), 161-70.

Note: This product is for in vitro research use only