

Eprazinone dihydrochloride

Catalog No: #S3034


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eprazinone dihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	453.44
Formulation	C ₂₄ H ₃₂ N ₂ O ₂ ·2HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
 SMILES: Cl.Cl.CCOC(CN1CCN(CC(C)C(=O)C2=CC=CC=C2)CC1)C1=CC=CC=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: Solubility: Water: 18 mg/mL (39.69 mM); DMSO: 1 mg/mL (2.2 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Thrall RS, et al. Exp Lung Res, 1992, 18(3), 409-420.

Note: This product is for in vitro research use only