

Cyclizine dihydrochloride

Catalog No: #S3037


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cyclizine dihydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Histamine Receptor antagonist
Calculated MW	339.3
Formulation	C18H24Cl2N2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: Cl.C1.CN1CCN(CC1)C(C1=CC=CC=C1)C1=CC=CC=C1
 Pathways: GPCR/G Protein;
 Immunology/Inflammation; Neuroscience
 Receptor: H1 receptor
 Boiling pt: Melting pt: Solubility: Water: 68 mg/mL (200.41 mM); DMSO: 7 mg/mL (20.63 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Vella-Brincat JW, et al. J Pain Symptom Manage. 2012 Mar;43(3):540-8.

Note: This product is for in vitro research use only