

Hederagenin

Catalog No: #S3039


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Hederagenin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	472.71
Formulation	C30H48O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@]12CC[C@@H]([C@@]([C@H]1CC[C@@]3[C@@H]2CC=C4[C@]3(CC[C@@]5([C@H]4CC(CC5)(C)C(=O)O)C)C)CO)O

Pathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lee CW, et al. Int Immunopharmacol. 2015 Dec;29(2):528-37.

Note: This product is for in vitro research use only