

Teriflunomide


Catalog No: #S3044

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Teriflunomide
Brief Description	Inhibitors
Purification	99.76%
Target Name	Dehydrogenase inhibitor
Calculated MW	270.21
Formulation	C12H9F3N2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: C/C(=C(\C#N)/C(=O)NC1=CC=C(C=C1)C(F)(F)F)/O
 Pathways: Metabolism
 Receptor: Dihydroorotate Dehydrogenase
 Boiling pt: Melting pt: Solubility: Ethanol: ,1 mg/mL ,1; Water:

References

1. Palmer AM. Curr Opin Investig Drugs. 2010 Nov;11(11):1313-23.

Note: This product is for in vitro research use only