

## Sancycline

Catalog No: #S3052


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sancycline
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Calculated MW	414.41
Formulation	C21H22N2O7
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OthersSMILES:CN(C)[C@H]1[C@@H]2C[C@@H]3CC4=C(C(=CC=C4)O)C(=C3C(=O)[C@@]2(C=C(C1=O)C(=O)N)O)OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mg/mLAppearance:PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Pinsuwan S, et al. Int J Pharm. 1999 Apr 20;181(1):31-40.;

Note: This product is for in vitro research use only