

Daurisoline

Catalog No: #S3054


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Daurisoline
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	610.75
Formulation	C37H42N2O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Other
 SMILES: CN1CCC2=CC(=C(C=C2[C@H]1CC3=CC=C(C=C3)OC4=C(C=CC(=C4)C[C@@H]5C6=CC(=C(C=C6CCN5C)OC)O)OC)O)OC
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: 96-102°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Liu QN, et al. Am J Chin Med. 2010;38(1):37-49.

Note: This product is for in vitro research use only