

Liensinine Perchlorate

Catalog No: #S3055


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Liensinine Perchlorate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	711.2
Formulation	C37H43ClN2O10
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: CN1CCC2=CC(=C(C=C2[C@H]1CC3=CC=C(C=C3)O)OC4=C(C=CC(=C4)C[C@@H]5C6=CC(=C(C=C6CCN5C)OC)OC)O)OC. OCl(=O)(=O)=O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Zhou J, et al. Autophagy. 2015;11(8):1259-79.

Note: This product is for in vitro research use only