

Bavachinin


Catalog No: #S3056

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bavachinin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	338.4
Formulation	C ₂₁ H ₂₂ O ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Other
SMILES: CC(=CCC1=C(C=C2C(=C1)C(=O)C[C@H](O2)C3=CC=C(C=C3)O)OC)C
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Chen X, et al. FEBS Lett. 2013 Sep 17;587(18):2930-5.

Note: This product is for in vitro research use only