


Description

Product Name	PD-173955
Brief Description	Inhibitors
Purification	99.28%
Target Name	Bcr-Abl inhibitor; Src inhibitor
Calculated MW	443.35
Formulation	C21H16Cl2N4OS
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: CSC1=CC=CC(NC2=NC=C3C=C(C(=O)N(C)C3=N2)C2=C(Cl)C=CC=C2Cl)=C1
 Pathways: Angiogenesis; Cytoskeletal Signaling
 Receptor: Bcr-Abl; Src
 Boling pt: Melting pt: Solubility: DMSO: 15 mg/mL warmed (33.83 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Wisniewski D, et al. Y Res. 2002, 62(15), 4244-4255.

Note: This product is for in vitro research use only