

GANT61

Catalog No: #S3070


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GANT61
Brief Description	Inhibitors
Purification	99.50%
Target Name	Hedgehog/Smoothened inhibitor
Calculated MW	429.6
Formulation	C27H35N5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:CancerSMILES:CN(C)C1=CC=CC=C1CN1CCCN(CC2=CC=CC=C2N(C)C)C1C1=CC=NC=C1Pathways:GPCR/G Protein; Stem CellsReceptor:GLI1Boiling pt:Melting pt:Solubility:Ethanol: 12 mg/mL (27.93 mM)Appearance:Faint Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Lauth M, et al. Proc Natl Acad Sci, 2007, 104(20), 8455-8560

Note: This product is for in vitro research use only