

## GLPG0634 analogue

Catalog No: #S3076


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	GLPG0634 analogue
Brief Description	Inhibitors
Purification	99.50%
Target Name	JAK inhibitor
Calculated MW	425.5
Formulation	C21H23N5O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cancer  
 SMILES: O=C(NC1=NN2C(C=CC=C2C=CC=C(CN3CCS(=O)(=O)CC3)C=C2)=N1)C1CC1  
 Pathways: Angiogenesis; Chromatin/Epigenetic; JAK/STAT signaling; Stem Cells  
 Receptor: JAK1; JAK2; JAK3  
 Boiling pt: Melting pt: Solubility: DMSO: 0.1 mg/mL (0.24 mM)  
 Appearance: Gray Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Van Rompaey L, et al. J Immunol. 2013, 191(7), 3568-3577.

Note: This product is for in vitro research use only