

Mesoridazine Besylate

Catalog No: #S3086


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mesoridazine Besylate
Brief Description	Inhibitors
Purification	99.27%
Target Name	5-HT Receptor antagonist; Dopamine Receptor antagonist
Calculated MW	544.75
Formulation	C27H32N2O4S3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous

system SMILES: OS(=O)(=O)C1=CC=CC=C1.CN1CCCC1CCN1C2=CC=CC=C2SC2=C1C=C(C=C2)S(C)=O Pathways: GPCR/G Protein;

Neuroscience Receptor: 5-HT; Dopamine Binding pt: Melting pt: Solubility: DMSO: 10 mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Im SH, et al. J Chromatogr B Analyt Technol Biomed Life Sci. 2014 May 1;958:117-23.

Note: This product is for in vitro research use only