

Tyramine


Catalog No: #S3095

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyramine
Brief Description	Inhibitors
Purification	98.00%
Target Name	MAO inhibitor
Calculated MW	137.18
Formulation	C ₈ H ₁₁ NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: NCCC1=CC=C(O)C=C1
Pathways: Metabolism
Receptor: MAOs
Boiling pt: Melting pt: Solubility: DMSO: 10 mM
Appearance: Yellow Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Grundemann D, et al. Nat Neurosci. 1998 Sep;1(5):349-51.

Note: This product is for in vitro research use only