

Potassium canrenoate

Catalog No: #S3096


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Potassium canrenoate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Glucocorticoid Receptor antagonist
Calculated MW	396.56
Formulation	C ₂₂ H ₂₉ KO ₄
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@]12CCC(=O)C=C1C=C[C@@H]3[C@@H]2CC[C@]4([C@H]3CC[C@]4(CCC(=O)[O-])O)C.[K+]

Pathways:Endocrinology/Hormones
 Receptor:Mineralocorticoid Receptor
 Boling pt:Melting pt:Solubility:DMSO: 10 mM
 Appearance:White Powder
 Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
 Stock solution can be stored below -20 °C for several month.

References

- Alexandre J, et al. J Cardiovasc Pharmacol Ther. 2015 May;20(3):313-21.

Note: This product is for in vitro research use only