

Mycophenolic Acid Impurity

Catalog No: #S3097


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mycophenolic Acid Impurity
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others inhibitor
Formulation	N/A
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others
SMILES: CC1=C(C=C(C2=C1COC2=O)O)O
Pathways: Others
Receptor: Others
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several month.

References

N/A

Note: This product is for in vitro research use only