


Description

Product Name	PRT 4165
Brief Description	Inhibitors
Purification	98.00%
Target Name	BMI-1 inhibitor; Others inhibitor
Calculated MW	235.24
Formulation	C15H9NO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Others
SMILES: C1=CC=C2C(=C1)C(=O)C(=CC3=CN=CC=C3)C2=O
Pathways: Others
Receptor: BMI-1; Ring1A
Boiling pt: Melting pt: Solubility: Soluble in DMSO
Appearance: Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ismail, et al. J Biol Chem. 2013 Sep 13;288(37):26944-54.

Note: This product is for in vitro research use only