

LY2835219 mesylate

Catalog No: #S3111


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	LY2835219 mesylate
Brief Description	Inhibitors
Purification	99.00%
Target Name	CDK inhibitor
Calculated MW	602.7
Formulation	C27H32F2N8CH4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer SMILES: CS(O)(=O)=O.CCN1CCN(CC2=CN=C(NC3=NC=C(F)C(=N3)C3=CC4=C(N=C(C)N4C(C)C)C(F)=C3)C=C2)CC1 Pathways: Cell Cycle/Checkpoint Receptor: CDK4; CDK6 Boiling pt: Melting pt: Solubility: Ethanol: 24 mg/mL (39.82 mM); Water: 100 mg/mL (165.92 mM); DMSO: 83 mg/mL (137.71 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

- Sanchez-Martinez et al. Mol Y Ther, 2011, 10(11 Suppl), Abstract nr B234.

Note: This product is for in vitro research use only