

## Verteporfin

Catalog No: #S3112


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Verteporfin
Brief Description	Inhibitors
Purification	99.00%
Target Name	VDA inhibitor
Calculated MW	718.79
Formulation	C41H42N4O8
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Cardiovascular

SystemSMILES:O=C(O)CCC1=C(C)C(/C=C2C(C=C)=C(C)/C(N/2)=C/3)=N/C1=C\C(N4)=C(CCC(OC)=O)C(C)=C4/C=C5[C@]6(C)[C@H](C(OC)=O)C(C(OC)=O)=CC=C6C3=N/5Pathways:AngiogenesisReceptor:VDABoling pt:Melting pt:Solubility:DMSO: 100 mg/mL (139.12 mM)Appearance:BlackPowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Schmidt-Erfurth U, et al. Surv Ophthalmol, 2000, 45(3), 195-214.

Note: This product is for in vitro research use only