

Foretinib (GSK1363089)

Catalog No: #S3113


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Foretinib (GSK1363089)
Brief Description	Inhibitors
Purification	99.80%
Target Name	c-Met/HGFR inhibitor; Tie-2 inhibitor; VEGFR inhibitor
Calculated MW	632.65
Formulation	C34H34F2N4O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Cancer
 SMILES: COC1=CC2=C(OC3=C(F)C=C(C=C3)N(C(=O)C3(CC3)C(N)=O)C3=CC=C(F)C=C3)C=CN=C2C=C1OCCCN1CCOCC1
 Pathways: Angiogenesis; Tyrosine Kinase/Adaptors Receptor: Met; RON; Tie-2; VEGFR2?(KDR); VEGFR3/FLT4
 Boiling pt: Melting pt: Solubility: DMSO: 127 mg/mL (200.74 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Qian F, et al. Y Res, 2009, 69(20), 8009-8016.

Note: This product is for in vitro research use only