

Povidone-iodine (PVP-I)

Catalog No: #S3114


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Povidone-iodine (PVP-I)
Brief Description	Inhibitors
Purification	I:9%-12%
Target Name	Others
Calculated MW	364.95
Formulation	C6H9I2NO
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:CCC(C)N1CCCC1=OPathways:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Brownness PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Akcay E, et al. Childs Nerv Syst. 2012, 28(12):2071-5.

Note: This product is for in vitro research use only