


Description

Product Name	JNJ-63533054
Brief Description	Inhibitors
Purification	98.00%
Target Name	GPR agonist
Calculated MW	316.78
Formulation	C17H17ClN2O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Others SMILES: CC(C)(c1ccccc1)NC(=O)NC(=O)c2ccc(Cl)cc2 Pathways: Endocrinology/Hormones Receptor: hGPR139 Boling pt: Melting pt: Solubility: Soluble in DMSO Appearance: Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Dvorak CA, et al. ACS Med Chem Lett. 2015 Jul 20;6(9):1015-8.

Note: This product is for in vitro research use only