

## Betulin

Catalog No: #S3121


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Betulin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Lipid inhibitor
Calculated MW	442.72
Formulation	C30H50O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Other SMILES: CC(=C)[C@H]1CC[C@]2([C@H]1[C@H]3CC[C@@H]4[C@]5(CC[C@@H](C([C@H]5CC[C@]4([C@@]3(CC2)C)C)(C)C)O)C)COP Pathways: Metabolism Receptor: SREBP (K562 cell) Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Zhao H, et al. Life Sci. 2016 Jan 1;144:185-93.

Note: This product is for in vitro research use only