

(+)-Fangchinoline

Catalog No: #S3122


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	(+)-Fangchinoline
Brief Description	Inhibitors
Purification	99.00%
Target Name	FAK inhibitor
Calculated MW	608.72
Formulation	C37H40N2O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Cardiovascular

SystemSMILES:CN1CCC2=CC(=C3C=C2[C@@H]1CC4=CC=C(C=C4)OC5=C(C=CC(=C5)C[C@H]6C7=C(O3)C(=C(C=C7CCN6C)OC)O)OC)OCPa

thways:Angiogenesis; Cytoskeletal SignalingReceptor:FAKBoling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:WhitePowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Guo B, et al. J Drug Target. 2015 Apr;23(3):266-74.

Note: This product is for in vitro research use only