

Bethanechol chloride

Catalog No: #S3126


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bethanechol chloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR agonist
Calculated MW	196.68
Formulation	C7H17N2O2Cl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous System
 SMILES: [Cl-].CC(C[N+](C)(C)C)OC(N)=O
 Pathways: Neuroscience
 Receptor: muscarinic receptor
 Boiling pt: Melting pt: Solubility: Ethanol: 39 mg/mL (198.29 mM); Water: 39 mg/mL (198.29 mM); DMSO: 39 mg/mL (198.29 mM)
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Harbaugh RE, et al. J Neurosurg. 1989 Oct;71(4):481-6.

Note: This product is for in vitro research use only