

## Glycinamide Hydrochloride

Catalog No: #S3127


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Glycinamide Hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	110.54
Formulation	C <sub>2</sub> H <sub>7</sub> CIN <sub>2</sub> O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: Cl.NCC(N)=O  
Pathways: Others  
Receptor: Others  
Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Beyer C, et al. Life Sci. 2013 Mar 21;92(10):576-81.

Note: This product is for in vitro research use only