

# Boldenone Undecylenate


Catalog No: #S3137

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Boldenone Undecylenate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Androgen Receptor agonist
Calculated MW	452.67
Formulation	C30H44O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

### Research

Area:CancerSMILES:C[C@]12CC[C@H]3[C@H]([C@@H]1CC[C@@H]2OC(=O)CCCCCCCC=C)CCC4=CC(=O)C=C[C@]34CPathways:Endocrinology/HormonesReceptor:Androgen ReceptorBoiling pt:Melting pt:Solubility:DMSO: 41 mg/mLAppearance:YellowLiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Oda SS, El-Ashmawy IM. Int J Exp Pathol. 2012 Jun;93(3):172-8.

Note: This product is for in vitro research use only