

## Astaxanthin

Catalog No: #S3138


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Astaxanthin
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	596.84
Formulation	C40H52O4
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

## Research

Area:OtherSMILES:CC1=C(C(C[C@@H](C1=O)O)(C)C)/C=C/C=C/C(=C/C=C/C(=C/C=C/C(=C/C=C/C(=O)[C@H](CC2(C)C)O)C)\C)/C/CPATHWAYS:OthersReceptor:OthersBoiling pt:Melting pt:Solubility:Soluble in ChloroformAppearance:RedPowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Zhang J, et al. J Cell Mol Med. 2015 Sep;19(9):2215-31.

Note: This product is for in vitro research use only