

Sultamicillin Tosylate

Catalog No: #S3139


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Sultamicillin Tosylate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Antibiotic inhibitor
Calculated MW	766.86
Formulation	C32H38N4O12S3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1=CC=C(C=C1)S(=O)(=O)O.CC1([C@@H](N2[C@H](S1)[C@@H](C2=O)NC(=O)[C@H](C3=CC=CC=C3)N)C(=O)OC(=O)[C@H]4C(S(=O)(=O)[C@H]5N4C(=O)C5)(C)C)C Pathways: Microbiology&virology Receptor: Antibiotic Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Friedel HA, et al. *Drugs*. 1989 Apr;37(4):491-522.

Note: This product is for in vitro research use only