

Bendazac L-lysine

Catalog No: #S3141


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bendazac L-lysine
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	428.48
Formulation	C22H28N4O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C1=CC=C(C=C1)CN2C3=CC=CC=C3C(=N2)OCC(=O)O.C(CCN)C[C@@H](C(=O)O)N Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Iuliano G, et al. Enzyme. 1989;42(4):235-7.

Note: This product is for in vitro research use only