

D-(+)-Trehalose dihydrate

Catalog No: #S3145


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	D-(+)-Trehalose dihydrate
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	378.33
Formulation	C ₁₂ H ₂₆ O ₁₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: C([C@@H]1[C@H]([C@@H]([C@H]([C@H](O1)O)[C@@H]2[C@@H]([C@H]([C@@H]([C@H](O2)CO)O)O)O)O)O.O.O Pathways: Others Receptor: Others Boiling pt: Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Sapir L, Harries D. J Phys Chem B. 2011 Feb 3;115(4):624-34.

Note: This product is for in vitro research use only