

PD1-PDL1 inhibitor 2

Catalog No: #S3146


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PD1-PDL1 inhibitor 2
Brief Description	Inhibitors
Purification	98.54%
Target Name	PD-1/PD-L1 inhibitor
Calculated MW	419.52
Formulation	C25H29N3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cancer
 SMILES: CC1=C(COC2=CC=C(CNCCNC(C)=O)C(OC)=N2)C=CC=C1C3=CC=CC=C3
 Pathways: Apoptosis; Cell Cycle/Checkpoint; Immunology/Inflammation
 Receptor: PD-1/PD-L1
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Louis S. Chupak, et al. Compounds useful as immunomodulators

Note: This product is for in vitro research use only